



Recombinant Human Splicing factor 1 (SF1)

Product Code	CSB-MP021120HU
Abbreviation	SF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15637
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ATGANATPL DFPSKKRKRS RWNQDTMEQK TVIPGMPTVI PPGLTREQER AYIVQLQIED LTRKLRTGDL GIPPNPEDRS PSPEPIYNSE GKRLNTREFR TRKKLEEERH NLITEMVALN PDFKPPADYK PPATRVSDKV MIPQDEYPEI NFVGLLIGPR GNTLKNIEKE CNAKIMIRGK GSVKEGKVGR KDGQMLPGED EPLHALVTAN TMENVKKAVE QIRNILKQGI ETPEDQNDLR KMQLRELARL NGTLREDDNR ILRPWQSSET RSITNTTVCT KCGGAGHIAS DCKFQRPGDP QSAQDKARMD KEYLSLMAEL GEAPVPASVG STSGPATTPL ASAPRPAAPA NNPPPPSLMS TTQSRPPWMN SGPSESRPYH GMHGGGGPGGP GGGPHSFPH LPSLTGGHGG HPMQHNPNGP PPPWMQPPPP PMNQGPHPG HHGPPMDQY LGSTPVGSGV YRLHQGKGM PPPPMGMMPP PPPPPSGQPP PPPSGPLPPW QQQQQQPPPP PPSSSMASST TPLPWQQNTT TTTTSAGTGS IPPWQQQQA AAASPGAPQM QGNPTMVPLP PGVQPPLPPG APPPPPPPP GSAGMMYAPP PPPPPPMDPS NFVTMMGMGV AGMPPFGMPP APPPPPPQN
Source	Mammalian cell
Target Names	SF1
Protein Names	Recommended name: Splicing factor 1 Alternative name(s): Mammalian branch point-binding protein Short name= BBP Short name= mBBP Transcription factor ZFM1 Zinc finger gene in MEN1 locus Zinc finger protein 162
Expression Region	2-639
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.