



Recombinant Human Zinc finger protein 174 (ZNF174)

Product Code	CSB-MP026555HU
Abbreviation	ZNF174
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15697
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAKMEITLS SNTEASSKQE RHIIAKLEEK RGPPLQKNCP DPCLCRQSFR RFCYQEVSGP QEALSQLRQL CRQWLQPELH TKEQILELLV MEQFLTILPP EIQARVRHRC PMSSKEIVTL VEDFHRASKK PKQWVAVCMQ GQKVLLEKTG SQLGEQELPD FQPQTPRRDL RESSPAEPSQ AGAYDRLSPH HWEKSPLLQE PTPKLAGTEA PRMRSDNKEN PQQEGAKGAK PCAVSAGRSK GNGLQNPEPR GANMSEPRLS RRQVSSPNAQ KPFAHYQRHC RVEYISSPLK SHPLRELKKS KGGKRSLSNR LQHLGHQPTR SAKKPYKCDD CGKSFTWNSE LKRHKRVHTG ERPYTCGECG NCFGRQSTLK LHQRIHTGEK PYQCGQCGKS FRQSSNLHQH HRLHHGD
Source	Mammalian cell
Target Names	ZNF174
Protein Names	Recommended name: Zinc finger protein 174 Alternative name(s): AW-1 Zinc finger and SCAN domain-containing protein 8
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.