



Recombinant Human Pancreatic lipase-related protein 3 (PNLIPRP3)

Product Code	CSB-BP018265HU
Abbreviation	PNLIPRP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17RR3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	KEV CYERLGCFKD GLPWTRTFST ELVGLPWSPE KINTRFLLYT IHNPNAVQEI SAVNSSTIQA SYFGTDKITR INIAGWKTDG KWQRDMCNVL LQLEDINCIN LDWINGSREY IHAVNNLRVV GAEVAYFIDV LMKKFEYSPS KVHLIGHSLG AHLAGEAGSR IPGLGRITGL DPAGPFFHNT PKEVRLDPSD ANFVDVIHTN AARILFELGV GTIDACGHLD FYPNGGKHP GCEDLITPLL KFNFNAYKKE MASFFDCNHA RSYQFYAESI LNPDAFIAYP CRSYTSFKAG NCFPCSKEGC PTMGHFADRF HFKNMKTNGS HYFLNTGSLS PFARWRHKLS VKLSGSEVTQ GTVFLRVGGA VRKTGEFAIV SGKLEPGMTY TKLIDADVNV GNITSVQFIW KKHLFEDSQN KLGAEMVINT SGKYGYKSTF CSQDIMGPNI LQNLKPC
Source	Baculovirus
Target Names	PNLIPRP3
Protein Names	Recommended name: Pancreatic lipase-related protein 3 Short name= PL-RP3 EC= 3.1.1.3
Expression Region	18-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.