



Recombinant Human Zinc finger protein 460 (ZNF460)

Product Code	CSB-BP026766HU
Abbreviation	ZNF460
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14592
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAWMAPAQ ESVTFEDVAV TFTQEEWGQL DVTQRALYVE VMLETGGLLV ALGDSTKPET VEPIPSHLAL PEEVSLQEQL AQGVPRYSYL GQAMDQDGPS EMQEYFLRPG TDPQSEKHLHG KMSLEHEGLA TADGICSM MI QNQVSPEDAL YGFDSYGPVT DSLIHEGENS YKFEEMFNEN CFLVQHEQIL PRVKPYDCPE CGKAFGKSKH LLQHHIHTG EKPYKLECG KDFNRRSHLT RHQRTHNGDK PFVCSECGRT FNRGSHLTRH QRVHSGEKPF VCNECGKRAFT YRSNFVLH NK SHNEKKPFAC SECGKGFYES TALIQHFIH TGERPFKCLE CGKAFNCRSH LKQHERIHTG EKPFVCSQCG KAFTHYSTYV LHERAHTGEK PFECKECKGA FSIRKDLIRH FNIHTGEKPY ECLQCGKRAFT RMSGLTRHQW IHTGEKPYVC IQCGKAF CRT TNLIRHFSIH TGEKPYECVE CGKAFNRRSP LTRHQRIHTA EKSHEPIQSG NVSCESTD LI QHSIIHTESS PVS AVN METP SIAAHSSSLD INGFIVEETL PL
Source	Baculovirus
Target Names	ZNF460
Protein Names	Recommended name: Zinc finger protein 460 Alternative name(s): Zinc finger protein 272 Zinc finger protein HZF8
Expression Region	1-562
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.