



Recombinant Human Nucleolar and coiled-body phosphoprotein 1 (NOLC1)

Product Code	CSB-BP015929HU
Abbreviation	NOLC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14978
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MADAGIRRVV PSDLYPLVLG FLRDNQLSEV ANKFAKATGA TQQDANASSL LDIYSFWLKS AKVPERKLQA NGPVAKKAKK KASSSDSEDS SEEEEEVQGP PAKKAAPPAK RVGLPPGKAA AKASESSSSE ESSDDDDDEED QKKQPVQKGV KPQAKAAKAP PPKAKSSDSD SDSSSEDEPP KNQKPKITPV TVKAQTKAPP KPARAAPKIA NGKAASSSSS SSSSSSSDSDS EEEKAAATPK KTVPKKQVVA KAPVKAATTP TRKSSSSSEDS SSDEEEEEQKK PMKNKPGPYS SVPPPSAPP KKSLGTQPPK KAVEKQQPVE SSEDSSDESD SSSEEEKKPP TKAVVSKATT KPPPAKKAEE SSSDSSDSDS SEDDEAPSKP AGTTKNSSNK PAVTTKSPAV KPAAAPKQPV GGGQKLLTRK ADSSSSEEEES SSSEEEKTKK MVATTKPKAT AKAALSLPAK QAPQGSRDSS SDSDSSSSSEE EEEKTSKSAV KKKPQKVAGG AAPSKPASAK KGKAESSNSS SSSDSSSEEEE EKLKGGKSPR PQAPKANGTS ALTAQNGKAA KNSEEEEEEEK KKAAVVVS KS GSLKRRKQNE AAKEAETPQA KKIKLQTPNT FPKRKKGEKR ASSPFRVRE EEIEVDSRVA DNSFDAKRG AGDWGERANQ VLKFTKGKSF RHEKTKKRG SYRGG SISVQ VNSIKFDSE
Source	Baculovirus
Target Names	NOLC1
Protein Names	Recommended name: Nucleolar and coiled-body phosphoprotein 1 Alternative name(s): 140 kDa nucleolar phosphoprotein Short name= Nopp140 Hepatitis C virus NS5A-transactivated protein 13 Short name= HCV NS5A-transactivated protein 1
Expression Region	1-699
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.