



# Recombinant Human Nucleolar and coiled-body phosphoprotein 1 (NOLC1)

<b>Product Code</b>	CSB-EP015929HU
<b>Abbreviation</b>	NOLC1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q14978
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MADAGIRRVV PSDLYPLVLG FLRDNQLSEV ANKFAKATGA TQQDANASSL LDIYSFWLKS AKVPERKLQA NGPVAKKAKK KASSSDSEDS SEEEEEVQGP PAKKAAPPAK RVGLPPGKAA AKASESSSSE ESSDDDEED QKKQPVQKGV KPQAKAAKAP PPKAKSSDSD SDSSSEDEPP KNQKPKITPV TVKAQTKAPP KPARAAPKIA NGKAASSSSS SSSSSSSDDS EEEKAAATPK KTVPKKQVVA KAPVKAATTP TRKSSSSEDS SSDEEEEQKK PMKNKPGPYS SVPPPSAPP KKSLGTQPPK KAVEKQQPVE SSEDSSDESD SSSEEEKKPP TKAVVSKATT KPPPAKAAE SSSDSSSDSDS SEDDEAPSKP AGTTKNSSNK PAVTTKSPAV KPAAAPKQPV GGGQKLLTRK ADSSSSEEEES SSSEEEKTKK MVATTKPKAT AKAALSLPAK QAPQGSRDSS SDSDSSSSEE EEEKTSKSAV KKKPQKVAGG AAPSKPASAK KGKAESSNSS SSDDSSSEEEE EKLKGGKSPR PQAPKANGTS ALTAQNGKAA KNSEEEEEEEK KKAAVVVS KS GSLKKRKQNE AAKEAETPQA KKIKLQTPNT FPKRKKGEKR ASSPFRVRE EEIEVDSRVA DNSFDAKRG AGDWGERANQ VLKFTKGKSF RHEKTKKRG SYRGG SISVQ VNSIKFDSE
<b>Source</b>	E.coli
<b>Target Names</b>	NOLC1
<b>Protein Names</b>	Recommended name: Nucleolar and coiled-body phosphoprotein 1 Alternative name(s): 140 kDa nucleolar phosphoprotein Short name= Nopp140 Hepatitis C virus NS5A-transactivated protein 13 Short name= HCV NS5A-transactivated protein 1
<b>Expression Region</b>	1-699
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.