



Recombinant Human Zinc finger protein 239 (ZNF239)

Product Code	CSB-MP026605HU
Abbreviation	ZNF239
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16600
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASTITGSQD CIVNHRGEVD GEPELDISPC QQWGEASSPI SRNRDSVMTL QSGCFENIES ETYLPLKVSS QIDTQDSSVK FCKNEPQDHQ ESRRLFVMEE STERKVIKGE SCSENQVVKL VSDGQELASP LLNGEATCQN GQLKESLDPI DCNCKDIHWG KSQVVSCSQQ RAHTEEKPCD HNNCGKILNT SPDGHPIYEKI HTAEKQYECS QCGKNFSQSS ELLHQRDHT EEKPYKCEQC GKGFTRSSSL LIHQAVHTDE KPYKCDKCGK GFTRSSSLI HHAVHTGEKP YKCDKCGKGF SQSSKLHIHQ RVHTGEKPYE CEECGMSFSQ RSNLHIHQRV HTGERPYKCG ECGKGFSQSS NLHIHRCIHT GEKPYQCYEC GKGFSQSSDL RIHLRVHTGE KPYHCGKCGK GFSQSSKLLI HQRVHTGEKP YECSKCGKGF SQSSNLHIHQ RVHKKDPR
Source	Mammalian cell
Target Names	ZNF239
Protein Names	Recommended name: Zinc finger protein 239 Alternative name(s): Zinc finger protein HOK-2 Zinc finger protein MOK-2
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.