



Recombinant Human Proline-rich nuclear receptor coactivator 1 (PNRC1)

Product Code	CSB-BP018289HU
Abbreviation	PNRC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12796
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTVVSVPQRE PLVLGGRLAP LGFSSRGYFG ALPMVTTAPP PLPRIPDPRA LPPTLFLPHF LGGDGPCLTP QPRAPAALPN RSLAVAGGTP RAAPKKRRKK KVRASPAGQL PSRFHQYQQH RPSLEGGRSP ATGPSGAQEV PGAAALAPS PAAAAGTEGA SPDLAPLRPA APGQTPLRKE VLKSKMGKSE KIALPHGQLV HGIHLYEQPK INRQKSKYNL PLTKITSKR NENFWQDSV SSDRIQKQEK KPFKNTENIK NSHLKSAFL TEVSQKENYA GAKFSDPPSP SVLPKPPSHW MGSTVENSQ NRELMVHLK TLLKVQT
Source	Baculovirus
Target Names	PNRC1
Protein Names	Recommended name: Proline-rich nuclear receptor coactivator 1 Alternative name(s): Proline-rich protein 2 Protein B4-2
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.