



# Recombinant Human Growth factor receptor-bound protein 10 (GRB10)

<b>Product Code</b>	CSB-BP009887HU
<b>Abbreviation</b>	GRB10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q13322
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MALAGCPDSF LHPYYQDKV EQTPRSQQDP AGPGLPAQSD RLANHQEDDV DLEALVNDMN ASLESLSYAC SMQSDTVPLL QNGQHARSQP RASGPPRSIQ PQVSPRQRVQ RSQPVHILAV RRLQEEDQQF RTSSLPAIPN PFPELCGPGS PPVLTGSLP PSQAAAKQDV KVFSEDGTSK VVEILADMTA RDLQQLLVYK SHCVDDNSWT LVEHHPHLGL ERCLEDHELV VQVESTMASE SKFLFRKNYA KYEFFKNPMN FFPEQMVTWC QQSNGSQTQL LQNFLNSSSC PEIQGFLHVK ELGKKS WKKL YVCLRRSGLY CSTKGTSKEP RHLQLLADLE DSNIFSLIAG RKQYNAPTDH GLCIKPNKVR NETKELRLLC AEDEQTRTCW MTAFRLLKYG MLLYQNYRIP QQRKALLSPF STPVRSVSEN SLVAMDFSGQ TGRVIENPAE AQSAALEEGH AWRKRSTRMN ILGSQSPLHP STLSTVIHRT QHWFHGRISR EESHRIKQQ GLVDGLFLLR DSQSNPKAFV LTLCHHQKIK NFQILPCEDD GQTFPSLDDG NTKFSDLIQL VDFYQLNKGV LPCKLKHHC RVAL
<b>Source</b>	Baculovirus
<b>Target Names</b>	GRB10
<b>Protein Names</b>	Recommended name: Growth factor receptor-bound protein 10 Alternative name(s): GRB10 adapter protein Insulin receptor-binding protein Grb-IR
<b>Expression Region</b>	1-594
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	The product of this gene belongs to a small family of adapter proteins that are known to interact with a number of receptor tyrosine kinases and signaling molecules. This gene encodes a growth factor receptor-binding protein that interacts with insulin receptors and insulin-like growth-factor receptors. Overexpression of some isoforms of the encoded protein inhibits tyrosine kinase



activity and results in growth suppression. This gene is imprinted in a highly isoform- and tissue-specific manner. Alternatively spliced transcript variants encoding different isoforms have been identified.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.