



# Recombinant Human Zinc finger protein 187 (ZNF187)

<b>Product Code</b>	CSB-EP026565HU
<b>Abbreviation</b>	ZNF187
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q16670
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MATALVSAHS LAPLNLKKEG LRVVREDHYS TWEQGFKLQG NSKGLGQEPL CKQFRQLRYE ETTGPREALS RLRELCQQWL QPETHTKEQI LELLVLEQFL IILPKELQAR VQEHHPESRE DVVVVLEDLQ LDLGETGQQD PDQPKKQKIL VEEMAPLKGV QEQQVRHECE VTKPEKEKGE ETRIENGKLI VVTDSCGRVE SSGKISEPME AHNEGSNLER HQAKPKEKIE YKCSEREQRF IQHLDLIEHA STHTGKKLCE SDVCQSSSLT GHKKVLSREK GHQCHECGKA FQRSSHLVRH QKIHLGKPY QCNECGKVFS QNAGLLEHLR IHTGEKPYLC IHCGKNFRRS SHLNRHQRIH SQEEPCECKE CGKTFSQALL LTHHQRIHSH SKSHQCNECG KAFSLTSDLI RHHRIHTGEK PFKCNICQKA FRLNSHLAQH VRIHNEEKPY QCSECGEAFR QRSGLFQHQR YHHKDKLA
<b>Source</b>	E.coli
<b>Target Names</b>	ZSCAN26
<b>Protein Names</b>	Recommended name: Zinc finger protein 187 Alternative name(s): Protein SRE-ZBP Zinc finger and SCAN domain-containing protein 26
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.