



Recombinant Human Zinc finger protein 187 (ZNF187)

Product Code	CSB-MP026565HU
Abbreviation	ZNF187
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16670
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MATALVSAHS LAPLNLKKEG LRVVREDHYS TWEQQGFKLQG NSKGLGQEPL CKQFRQLRYE ETTGPREALS RLRELCQQWL QPETHTKEQI LELLVLEQFL IILPKELQAR VQEHHPESRE DVVVVLEDLQ LDLGETGQQD PDQPKKQKIL VEEMAPLKGV QEQQVRHECE VTKPEKEKGE ETRIENGKLI VVTDSCGRVE SSGKISEPME AHNEGSNLER HQAKPKEKIE YKCSEREQRF IQHLDLIEHA STHTGKKLCE SDVCQSSSLT GHKKVLSREK GHQCHECGKA FQRSSHLVRH QKIHLGKPY QCNECGKVFS QNAGLLEHLR IHTGEKPYLC IHCGKNFRRS SHLNRHQRIH SQEEPCECKE CGKTFSQALL LTHHQRIHSH SKSHQCNECG KAFSLTSDLI RHHRIHTGEK PFKCNICQKA FRLNSHLAQH VRIHNEEKPY QCSECGEAFR QRSGLFQHQR YHHKDKLA
Source	Mammalian cell
Target Names	ZSCAN26
Protein Names	Recommended name: Zinc finger protein 187 Alternative name(s): Protein SRE-ZBP Zinc finger and SCAN domain-containing protein 26
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.