



# Recombinant Human Heterogeneous nuclear ribonucleoprotein A0 (HNRNPA0)

<b>Product Code</b>	CSB-EP010599HU-B
<b>Abbreviation</b>	HNRNPA0
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q13151
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MENSQLCKLF IGGLNVQTSE SGLRGHFEAF GTLTDCVVVV NPQTKRSRCF GFVTYSNVEE ADAAMAASPH AVDGNTVELK RAVSREDSAR PGAHAKVKKL FVGGLKGDVA EGDIEHFSQ FGTVEKAEII ADKQSGKKRG FGFVYFQNHD AADKAAVVKF HPIQGHRVEV KKAVPKEDIY SGGGGGGSRS SRGGRGGRGR GGGRDQNGLS KGGGGGYNSY GYGGGGGGGG YNAYGGGGGG SSGGSDYGN GFGGFGSYSQ HQSSYGPMKS GGGGGGGGSS WGGRSNSGPY RGGYGGGGGY GGSSF
<b>Source</b>	E.coli
<b>Target Names</b>	HNRNPA0
<b>Protein Names</b>	Recommended name: Heterogeneous nuclear ribonucleoprotein A0 Short name= hnRNP A0
<b>Expression Region</b>	1-305
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene belongs to the A/B subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. This protein has two repeats of quasi-RRM domains that bind RNAs, followed by a glycine-rich C-terminus.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.