



Recombinant Human Ras and Rab interactor 1 (RIN1)

Product Code	CSB-MP019726HU
Abbreviation	RIN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13671
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MESPGESGAG SPGAPSPSSF TTGHLAREKP AQDPLYDVPN ASGGQAGGPQ RPGRVVSRLRE RLLLTRPVWL QLQANAAAAL HMLRTEPPGT FLVRKSNTRQ CQALCMRLPE ASGPSFVSSH YILESPGGVS LEGSELMFPD LVQLICAYCH TRDILLPLQ LPRAIHHAAT HKELEAISHL GIEFWSSSLN IKAQRGPAGG PVLPQLKARS PQELDQGTGA ALCFFNPLFP GDLGPTKREK FKRSFKVRVS TETSSPLSPP AVPPPPVPVL PGAVPSQTER LPPCQLLRRE SSVGYRVPAG SGPSLPPMPS LQEVDGCGSPS SSEEELVPGS RGPATSPHL GRRRPLLRSM SAAFCSLLAP ERQVGRAAAA LMQDRHTAAG QLVQDLLTQV RAGPEPQELQ GIRQALSRAR AMLSAELGPE KLLSPKRLEH VLEKSLHCSV LKPLRPILAA RLRRRLAADG SLGRLAELR LARAQGPAGF GSHLSLPSPV ELEQVRQKLL QLLRTYSPSA QVKRLLQACK LLYMALRTQE GEGAGADEFL PLLSLVLAHC DLPELLLEAE YMSELLEPSL LTGEGGYLT SLSASLALLS GLGQAHTLPL SPVQELRRSL SLWEQRRLPA THCFQHLLRV AYQDPSSGCT SKTLAVPPEA SIATLNQLCA TKFRVTQPNT FGLFLYKEQG YHRLPPGALA HRLPTTGYLV YRRAEWPETQ GAVTEEEGSG QSEARSRGEE QGCQGDGDAG VKASPRDIRE QSETTAEGGQ GQAQEGPAQP GEPEAEGSRA AEE
Source	Mammalian cell
Target Names	RIN1
Protein Names	Recommended name: Ras and Rab interactor 1 Alternative name(s): Ras inhibitor JC99 Ras interaction/interference protein 1
Expression Region	1-783
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.