



Recombinant Human Zinc finger protein 192 (ZNF192)

Product Code	CSB-EP026568HU-B
Abbreviation	ZNF192
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15776
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEESRKPSAPSPPDQTPEEDLVIVKVEEDHGWDQESSLHESNPLGQEVFRL RFRQLRYQ ETLGPREALIQLRALCHQWLRPDLNTEQILELLVLEQFLTILPEELQTLVKEHQ LENGE EVTLLLEDLERQIDILGRPVSARVHGHRVLWEEVVHSASAPEPPNTQLQSEAT QHKSPVP QESQERAMSTSQSPTRSQKGSSGDQEMTATLLTAGFQTLEKIEDMAVSLIREE WLLDPSQ KDLCDNRPENFRNMFSLGGETRSENRELASKQVISTGIQPHGETAAKCNGD VIRGLEHE EARDLLGRLERQRGNPTQERRHKCDECGKSFAQSSGLVRHWRIHTGEKPYQ CNVCGKAFS YRSALLSHQDIHNKVKRYHCCEGKAFSQNTGLILHQRIHTGEKPYQCNQCGK AFSQSAG LILHQRIHSGERPYPECNECGKAFSSHSLIGHQRIHTGEKPYECDECGKTFRR SSHIGH QRSHTGEKPYKCNECGRAFSQKSLIEHQRIHTGERPYKCKEKGKAFNGNTG LIQHLRIH TGEKPYQCNECGKAFIQRSSLRHQRIHSGEKSEISV
Source	E.coli
Target Names	ZKSCAN8
Protein Names	Recommended name: Zinc finger protein 192 Alternative name(s): LD5-1 Zinc finger protein with KRAB and SCAN domains 8
Expression Region	1-578
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.