



Recombinant Bovine Transcription factor Dp-1 (TFDP1)

Product Code	CSB-BP023424BO
Abbreviation	TFDP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17QZ4
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	AKDAGLIEA NGELKVFIDQ NLSPGKGVVS LVAVHPSTVN SLGKQLLPKT FGQSNVNITQ QVVIGTPQRP AAPNTIVVGS PHTPNTHFVS QNQPSPDPSPW SAGKRNKGE KNGKGLRHFS MKVCEKVQRK GTTSYNEVAD ELVAEFSAAD SHILPSESAY DQKNIRRRVY DALNVLMMAMN IISKEKKEIK WIGLPTNSAQ ECQSLEVERQ RRLERIKQKQ SQLQELILQQ IAFKNLVQRN RQVEQQASRP PPNSVIHLP FIIVNTSKKT VIDCSISNDK FEYLFNFDNT FEIHDDIEVL KRMGMACGLE SGSCSPEDLR VARSLVPKAL EPYVTEMAQG SLGGVVFVASA VSTSNGTRLS ASDLANGADG ALATSSSGSQ YSGSRVETPV SCVGEDDEDD EDFNENEED
Source	Baculovirus
Target Names	TFDP1
Protein Names	Recommended name: Transcription factor Dp-1
Expression Region	2-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.