



Recombinant Streptococcus pyogenes serotype M4 Uracil-DNA glycosylase (ung)

Product Code	CSB-MP025641SAAW
Abbreviation	ung
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1J724
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M4 (strain MGAS10750)
Purity	>85% (SDS-PAGE)
Sequence	MTHSIWHEKI KSFLPEHYYG RINHFLDEAY ASGLVYPQRE NVFKALQVTP LEETKVLILG QDPYHGPKQA QGLSFSVPEE ISAPPSLINI LKELADDIGP RDHHDSTWA SQGVLLLNAC LTVPAGQANG HAGLIWEPFT DAVIKVLNEK DSPVVFILWG AYARKKKAFI TNSKHIIIES PHPSPLSSYR GFFGSKPFSR TNAILEKEGM TGVDWLK
Source	Mammalian cell
Target Names	ung
Protein Names	Recommended name: Uracil-DNA glycosylase Short name= UDG EC= 3.2.2.27
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.