



# Recombinant Adenylosuccinate lyase (R06C7.5)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP001397CXY                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Abbreviation</b>      | R06C7.5                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                        |
| <b>Uniprot No.</b>       | Q21774                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Immunogen Species</b> | Caenorhabditis elegans                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Sequence</b>          | MASEDKFESV LSTRYCKNSP LVSILSETNK ATLWRQLWIW LAEAEKELGL<br>KQVTQDAIDE MKSNRDVFDW PFIRSEERKL KHDVMAHNHA FGKLCPTAAG<br>IIHLGATSCF VQDNADLIAY RDSIDHILKR FATVIDRLAA FSLKNKEVVT<br>VGRTHYQTAS LVTVGKRGVL WAQELLMAFQ SLSEFRDKMR FRGIKGATGT<br>QDSFLTTFAG DESKVEALDE LVTKKANFSN RFLITGQTYS RQQDSQLVFS<br>LSLLGAAAKK VCTDIRVLQA FGELLEPF EK DQIGSSAMPY KKNPMKSERC<br>CALSRKLINA PQEALTILAD QGLERTLDDS AGRRMLIPDV LLTAEALLTT<br>LQNIFEGLSV QTDNVKKIVE DEIAFLGLEK AMMMLTEEGV DRQQAHAVIR<br>KTALEAKQLQ ATQKVDIRQT MADPFFDSVR DRVVGLVNNP INFTGRCVSQ<br>TESFIAKELK PTIDKYLDKS AGNVQLDV |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Target Names</b>      | adsl-1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Protein Names</b>     | Recommended name: Adenylosuccinate lyase Short name= ASL EC= 4.3.2.2<br>Alternative name(s): Adenylosuccinase Short name= ASase                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Expression Region</b> | 1-478                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                               |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                        |