



# Recombinant Eukaryotic translation initiation factor 5 (C37C3.2)

<b>Product Code</b>	CSB-EP007572CXY-B
<b>Abbreviation</b>	C37C3.2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q22918
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MALNVNRAVA DPFYRYKMPK LSAKVEGKGN GIKTVISNMS EIAKALERPP MYPTKYFGCE LGAQTNFDAK NERYIVNGEH DANKLQDILD GFIKKFVLCK SCENPETQLF VRKNNIKSKC KACGCSFDID LKHKLSTFIM KNPPKIDVDF SKAEQKNGKK TSGADAAA AV AADIIHNSDK GSSNDDDDDD WEPEPVEPNG MLSAGMGKLV LDKDLEKSEE QRLDMLHTFF LKAKEEDRIS DAKGQTALRD EAERLELKQK ASLLLLANVFL DEKVITDKQI SKHRNLLLRFLNDKKAQRY LLGGVEQVIH KHEAELLSKS AHIKSLYDE DVCEEDSLIS WGEKPSSKYV SKSFAKKIIE NSQPVLNWLK EAEETEES DDEIAFGGDV KESEFLRQQK EKAAREAQQK SAKATNGNAA AASGANDEED LDIDDI
<b>Source</b>	E.coli
<b>Target Names</b>	C37C3.2
<b>Protein Names</b>	Recommended name: Eukaryotic translation initiation factor 5 Short name= eIF-5
<b>Expression Region</b>	1-436
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.