



Recombinant Eunectes notaeus Brain-derived neurotrophic factor (BDNF)

Product Code	CSB-MP002655EPK
Abbreviation	BDNF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1X709
Product Type	Recombinant Protein
Immunogen Species	Eunectes notaeus (Yellow anaconda)
Purity	≥85% (SDS-PAGE)
Sequence	HSDPARRGEL SVCDSTSEWV TAAEKKTAVD MSGATVTVLE KVPVPKGQLK QYFYETKCSS KGYAKEGCRG IDKRYWNSQC RTTQSFVRAL TMDNKKRVG
Source	Mammalian cell
Target Names	BDNF
Protein Names	Recommended name: Brain-derived neurotrophic factor Short name= BDNF
Expression Region	101-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.