



Recombinant Dog Fibroblast growth factor 5 (FGF5)

Product Code	CSB-MP008632DO
Abbreviation	FGF5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20FD0
Product Type	Recombinant Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Purity	≥85% (SDS-PAGE)
Sequence	HGEKHLAPKG QPGPAATGRN PGGAGGSSTS GGTSSSSSSS VSSAPGASPG IRGSGSEQGS FQWSPSGRRT GSLYCRVGIG FHLQIYPDGK VNGSHEANML SILEIFAVSQ GIVGIRGVFS NKFLAMSKKG KLHASAKFTD DCKFRERFQE NSYNTYASAI HRSEPAGREW YVALNKRKGA KRGCSPRVKP QHVSTHFLPR FKQLEHPELS FTVTVPEKKK PPSHVKPKVP LSAPRKSPNT VKYRLKFRFG
Source	Mammalian cell
Target Names	FGF5
Protein Names	Recommended name: Fibroblast growth factor 5 Short name= FGF-5
Expression Region	21-270
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This gene was identified as an oncogene, which confers transforming potential when transfected into mammalian cells. Targeted disruption of the homolog of this gene in mouse resulted in the phenotype of abnormally long hair, which suggested a function as an inhibitor of hair elongation. Alternatively spliced transcript variants encoding different isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.