



Recombinant Papiro hamadryas Transcription factor 12 (TCF12)

Product Code	CSB-MP023293ERJ
Abbreviation	TCF12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q28772
Product Type	Recombinant Protein
Immunogen Species	Papiro hamadryas (Hamadryas baboon)
Purity	≥85% (SDS-PAGE)
Sequence	NKEKDENLHE PPSSDDMKSD DESSQKDIKV SSRGRTSTNE DEDLNPEQKI EREKERRMAN NARERLRVRD INEAFKELGR MCQLHLKSEK PQTLLILHQ AVAVILSLEQ QVRERLNPK AACLKRREEE KVSVSAEPP TTLPGTHPGL SETTNPNGHM
Source	Mammalian cell
Target Names	TCF12
Protein Names	Recommended name: Transcription factor 12 Short name= TCF-12 Alternative name(s): DNA-binding protein HTF4 E-box-binding protein Transcription factor HTF-4
Expression Region	1-160
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.