



# Recombinant Escherichia coli Argininosuccinate lyase (argH)

<b>Product Code</b>	CSB-EP002213EGW-B
<b>Abbreviation</b>	argH
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1R3V1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain UTI89 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MALWGGRFTQ AADQRFKQFN DSLRFDYRLA EQDIVGSVAW SKALVTVGVL TAEFAQLEE ALNVLLEDVR ARPQQILESD AEDIHSWVEG KLIDKVGQLG KKLHTGRSRN DQVATDLKLW CKDTVSELLT ANRQLQSALV ETAQNNQDAV MPGYTHLQRA QPVTFAHWCL AYVEMLARDE SRLQDALKRL DVSPGCGAL AGTAYEIDRE QLAGWLG FAS ATRNSLDSVS DRDHVLELLS AAAIGMVHLS RFAEDLIFFN TGEAGFVELS DRVTSGSSLM PQKKNPDAL LIRGKCGRVQ GALTGMMMTL KGLPLAYNKD MQEDKEGLFD ALDTWLDCLH MAALVLDGIQ VKRPRCQEA QQGYANATEL ADYLVAKGVP FREAHHIVGE AVVEAIRQ GK ALEDLPLSEL QKFSQVIGED VYPILSLQSC LDKRAAKGGV SPQQVAQAIA FAQARLG
<b>Source</b>	E.coli
<b>Target Names</b>	argH
<b>Protein Names</b>	Recommended name: Argininosuccinate lyase Short name= ASAL EC= 4.3.2.1 Alternative name(s): Arginosuccinase
<b>Expression Region</b>	1-457
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.