



Recombinant Ursus thibetanus Leptin (LEP)

Product Code	CSB-EP012870UAY
Abbreviation	LEP
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q1XG29
Product Type	Recombinant Protein
Immunogen Species	Ursus thibetanus (Asiatic black bear) (Selenarctos thibetanus)
Purity	>85% (SDS-PAGE)
Sequence	VPIRKVQDD TKTLIKTIVT RINDISHTQA VSSKQRVAGL DFIPGLHPVL SLSRMDQTLA IYQQILTSLH SRNVVQISND LENLRDLLHL LASSKSCPLP RARGLESFES LGGVLEASLY STEVVALSRL QAALQDMLRR LDLSPGC
Source	E.coli
Target Names	LEP
Protein Names	Recommended name: Leptin Alternative name(s): Obesity factor
Expression Region	22-167
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>