



Recombinant Palmitoyl-protein thioesterase 1 (ppt-1)

Product Code	CSB-EP018587CXY
Abbreviation	ppt-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20390
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	EFRN ATKQVPVVMW HGMGDCCCNP LSMGSKKLF EEQIPGVYVH SLQLGSSITK DIEHGFYANT NELVYMACIK IKNDPELKNG YNAIGFSQGA QFLRAVAQRC PNPPMKNLVS VGGQHQQVFG APYICIGNIM CNGVRRIDL GAYLPFVQKR VVQAQYWHDP NQVEEYKKRS IFLADINNEN NNNPTYKRNL LSLKNLVLVK FNQDHMVVVK DSSWFGFYKD GDIDTILPMN ETDLYKEDRI GLKKLHESGR IHFMDVDGDH LQIPRSVLVN DIIKKYFM
Source	E.coli
Target Names	ppt-1
Protein Names	Recommended name: Palmitoyl-protein thioesterase 1 Short name= PPT-1 EC= 3.1.2.22 Alternative name(s): Palmitoyl-protein hydrolase 1
Expression Region	17-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.