



# Recombinant Pig Clusterin (CLU)

<b>Product Code</b>	CSB-MP005595PI
<b>Abbreviation</b>	CLU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q29549
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	KE LQEMSTEGSK YVNKEIKNAL KEVKQIKTLI EQSNEERKSL LSSLEEAKKK KEDALNDTRD TETKLKGSQG LCNETMMALW EECKPCLKQT CMKFYARVCR SGSGLVGHQL EEFLNQSSPF YFWINGDRID SLMENDRQQS HVMDIMEDSF NRASNIMDEL FQDRFFNREP FDTQFFSPFG SSHRGSLLFFN PKSRFARNIM PFPLFTDLNY HDMFQPFFDM IHQAQQAMDA HLHRIPYHFP EAGVPENSND RAVCKEIRHN STGCLRMKDQ CEKCREILSV DCSASNSSQM QLRQELYTSL QMAEKFSKLY DQLLQSYQQK MLNTSSLLKQ LNEQFSWVSQ LANLTQNDDR YYLQVTTVNS HGSDPSVPSG LTKVVVKLFD SYPITLIIPQ EVSDPKFMET VAEEALQQYR QRSREE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CLU
<b>Protein Names</b>	Recommended name: Clusterin Alternative name(s): CP40 Complement cytolysis inhibitor Short name= CLI Cleaved into the following 2 chains: 1. Clusterin beta chain 2. Clusterin alpha chain
<b>Expression Region</b>	29-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.