



Recombinant Mouse Speckle-type POZ protein-like (Spopl)

Product Code	CSB-EP022602MO-B
Abbreviation	Spopl
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2M2N2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSREPTPLPGDMSTSPVAESWCYTQVKVVKFSYMWTINNFSFCREEMGEV LKSSTFSSG PNDKMKWCLRVNPKGLDDESKDYLSLYLLLVSCPKEVRAKFKFSLLNDKREE TKAMESQ RAYRFVQGGKDWGFKKFIRRDPELLDEANGLLPDDKLTLCFCEVSVVQDSVNVSG HTSTNTLK VPECRLAEDLGNLWENTRFTDCCFFVRGKEFKAHKSVLAARSPVFNAMFEHE MEECTKNR VEINDLDPEVFKEMMRFVYTGKAPNLDMADNLLAAADKYALERLKVMECEAL CSNLSVE NVADTLVLADLHSAEQLKAQAIDFINRCSVLRQLGCKDGNWNNNQATDIMET SGWKSMI QSHPHLVAEAFRALASSQCPQFGIPRKRLKQS
Source	E.coli
Target Names	Spopl
Protein Names	Recommended name: Speckle-type POZ protein-like Alternative name(s): HIB homolog 2
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.