



Recombinant Horse Clusterin (CLU)

Product Code	CSB-MP005595HO
Abbreviation	CLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q29482
Product Type	Recombinant Protein
Immunogen Species	Equus caballus (Horse)
Purity	>85% (SDS-PAGE)
Sequence	DKAVSDREL QEMSTQGSNY INKEIKNALK GVKQIKNLIE QTNEERKSL GTLEEAKKKK EGALNDTKDS EMKLIKESQGV CNETMTALWE ECKPCLKQTC MKFYARVCRS GSGLVGHQLE EFLNQSSPFY FWINGDRIDS LLENDRQQTH VLDVMQDSFD RASSIMDELF QDRFFTREPPQ DTYYYSPFSS PHRRSSLLFN PKSRFARNIM HFPMYRHLNF NDMFQPFDDM IHQAQQAMNL HLHRLPDQLP MTEFSEGDNH DRTVCKEIRH NSTGCLKMKD QCEKCQEILS VDCSTNNPSQ MQLRQELNNS LQLAEKFTKL YDELLQSYQE KMLNTSSLLK QLNEQFSWVS QLANLTQGED QYYLQVTTVS SHNSDSEVPS GLTRVVVKLF DSYPIVTVV EVVSRNPKF METVAEKALQ EYRQKNREE
Source	Mammalian cell
Target Names	CLU
Protein Names	Recommended name: Clusterin Cleaved into the following 2 chains: 1. Clusterin beta chain 2. Clusterin alpha chain
Expression Region	22-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.