



Recombinant *Macropus giganteus* Transthyretin (TTR)

Product Code	CSB-MP025270MPF
Abbreviation	TTR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q29616
Product Type	Recombinant Protein
Immunogen Species	<i>Macropus giganteus</i> (Eastern gray kangaroo)
Purity	>85% (SDS-PAGE)
Sequence	VHHESEHSKCPLMVKVLDAVRGRPAVNVDVKVFKKTEEQTWELFAAGKTNDN GEIHELTT DDKFGELYKVEFDTISYWKALGVSPFHEYADVFTANDAGHRHYTIAAQLSP YSFSTTA IVSNPTE
Source	Mammalian cell
Target Names	TTR
Protein Names	Recommended name: Transthyretin Alternative name(s): Prealbumin
Expression Region	23-149
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.