



# Recombinant Human Zinc finger protein 800 (ZNF800)

<b>Product Code</b>	CSB-BP027053HU
<b>Abbreviation</b>	ZNF800
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2TB10
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPLRDKYCQT DHHHHGCCEP VYILEPGDPP LLQQPLQTSK SGIQQIIECF RSGTKQLKHI LLKDVDTIFE CKLCRSLFRG LPNLITHKKF YCPPSLQMDD NLPDVNDKQS QAINDLLEAI YPSVDKREYI IKLEPIETNQ NAVFYISRT DNPIEVTSS STPEQTEVQI QETSTEQSKT VPVTDTEVET VEPPPVEIVT DEVAPTSDEQ PQESQADLET SDNSDFGHQL ICCLCRKEFN SRRGVRRHIR KVHKKMEEL KKYIETRKNP NQSSKGRSKN VLVPLSRSCP VCKKSFATKA NVRRHFDEVH RGLRRDSITP DIATKPGQPL FLDSISPKKS FKTRKQKSSS KAEYNLTACK CLLCKRKYSS QIMLKRHMQI VHKITLSGTN SKREKGPNT ANSSEIKVKV EPADSVESP PSITHSPQNE LKGTNHSNEK KNTPAAQKNK VKQDSESPKS TSPSAAGGQQ KTRKPKLSAG FDFKQLYCKL CKRQFTSKQN LTKHIELHTD GNNIYVKFYK CPLCTYETRR KRDVIRHITV VHKKSSRYLG KITASLEIRA IKKPIDFVLN KVAKRGPSRD EAKHSDSKHD GTSNSPSKKY EVADVGIEVK VTKNFSLHRC NKCGKAFACK TYLEHHKKTH KANASNSPEG NKTKGRSTRS KALV
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF800
<b>Protein Names</b>	Recommended name: Zinc finger protein 800
<b>Expression Region</b>	1-664
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.