



Recombinant Human Zinc finger protein 800 (ZNF800)

Product Code	CSB-EP027053HU
Abbreviation	ZNF800
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2TB10
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPLRDKYCQT DHHHHGCCEP VYILEPGDPP LLQQPLQTSK SGIQQIIECF RSGTKQLKHI LLKDVDTIFE CKLCRSLFRG LPNLITHKKF YCPPSLQMDD NLPDVNDKQS QAINDLLEAI YPSVDKREYI IKLEPIETNQ NAVFYISRT DNPIEVTESS STPEQTEVQI QETSTEQSKT VPVTDTEVET VEPPPVEIVT DEVAPTSDEQ PQESQADLET SDNSDFGHQL ICCLCRKEFN SRRGVRRHIR KVHKKMEEL KKYIETRKNP NQSSKGRSKN VLVPLSRSCP VCKKSFATKA NVRRHDFEVH RGLRRDSITP DIATKPGQPL FLDSISPKKS FKTRKQKSSS KAEYNLTACK CLLCKRKYSS QIMLKRHMQI VHKITLSGTN SKREKGPNT ANSSEIKVKV EPADSVESP PSITHSPQNE LKGTNHSNEK KNTPAAQKNK VKQDSESPKS TSPSAAGGQQ KTRKPKLSAG FDFKQLYCKL CKRQFTSKQN LTKHIELHTD GNNIYVKFYK CPLCTYETRR KRDVIRHITV VHKKSSRYLG KITASLEIRA IKKPIDFVLN KVAKRGPSRD EAKHSDSKHD GTSNSPSKKY EVADVGIEVK VTKNFSLHRC NKCGKAFACK TYLEHHKKTH KANASNSPEG NKTGGRSTRS KALV
Source	E.coli
Target Names	ZNF800
Protein Names	Recommended name: Zinc finger protein 800
Expression Region	1-664
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.