



# Recombinant Arabidopsis thaliana Mitogen-activated protein kinase 7 (MPK7)

<b>Product Code</b>	CSB-BP013465DOA
<b>Abbreviation</b>	MPK7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q39027
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAMLVEPPNG IKQQGKHYY S MWQTLFEIDT KYVPIKPIGR GAYGVVCS S I NRETNERVAI KKIHNVFENR VDALRTLREL KLLRHVRHEN VIALKDVMLP ANRSSFKDVY LUYELMDTDL HQIIKSSQSL SDDHCKYFLF QLLRGLKYLH SANILHRDLK PGNLLVNANC DLKICDFGLA RTSQGNEQFM TEYV VTRWYR APELLLCCDN YGTSIDVWSV GCIFAEILGR KPIFPGTECL NQLKLIINVV GSQQESDIRF IDNPKARRFI KSLPYSRGTH LSNLYPQANP LAIDLLQRML VFDPTKRISV TDALLHPYMA GLFDPGSNPP AHVPISLDID ENMEEPVIRE MMWNEMLYYH PEAEISNA
<b>Source</b>	Baculovirus
<b>Target Names</b>	MPK7
<b>Protein Names</b>	Recommended name: Mitogen-activated protein kinase 7 Short name= AtMPK7 Short name= MAP kinase 7 EC= 2.7.11.24
<b>Expression Region</b>	1-368
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.