



Recombinant Human Zinc finger protein 699 (ZNF699)

Product Code	CSB-EP026960HU-B
Abbreviation	ZNF699
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q32M78
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEERKTAEL QKNRIQDSVV FEDVAVDFTQ EEWALLDLAQ RNLYRDVMLE NFQNLASLGY PLHTPHLISQ WEQEEDLQTV KRELIQGIFM GEHREGFETQ LKTNESVASQ DICGEKISNE QKIVRFKRND SWFSSLHENQ ESCGIDYQNK SHERHLRNHM VENIYECYEE NQDGQTFSQV PNLDLKRNT EVKSCECHEC GKAFVDHSSL KSHIRSHTGS KPYQCKECKG AFHFLACFKK HMKTPTEEKP YECKECTKAF SCSSFFRAHM KIHIGKTYE CKECGKGFSC SSSLTEHKRI HSGDKPYECK ECGKAFSCSS SLSKHKRIHS GDKPYECKEC GKAFSSSSHL IIHIRIHTGE KPYECKECKG AFSESSKLTV HGRTHTGEKP YKCKECKGAY NCPSSLSIHM RKHTGEKPYE CLECGKAFYL PTSLNTHVKN QSREKPYECK ECGKAFSCPS SFRAHVRDHT GKIQYECKEC GKTFSRSSSL TEHLRTHSGE KPYECKECKG AFISSSHLTV HIRHTHTGEKP YECKKCGKAF IYPSALRIHM RTHHTGEKPYE CKECGKAFRH SSYLT VHARM HTGEKPFCECL ECGKAFSCPS SFRRHVRSHHT GEKPYECKEC GKAFVCPAYF RRHVKTHTRE NI
Source	E.coli
Target Names	ZNF699
Protein Names	Recommended name: Zinc finger protein 699 Alternative name(s): Hangover homolog
Expression Region	1-642
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.