



# Recombinant Human Zinc finger protein 699 (ZNF699)

<b>Product Code</b>	CSB-MP026960HU
<b>Abbreviation</b>	ZNF699
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q32M78
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEEERKTAEL QKNRIQDSVV FEDVAVDFTQ EEWALLDLAQ RNLYRDVMLE NFQNLASLGY PLHTPHLISQ WEQEEDLQTV KRELIQGIFM GEHREGFETQ LKTNESVASQ DICGEKISNE QKIVRFKRND SWFSSLHENQ ESCGIDYQNK SHERHLRNHM VENIYECYEE NQDGQTFSQV PNLDLKRNT EVKSCECHEC GKAFVDHSSL KSHIRSHTGS KPYQCKECKG AFHFLACFKK HMKTPTEEKP YECKECTKAF SCSSFFRAHM KIHIGKTYE CKECGKGFSC SSSLTEHKRI HSGDKPYECK ECGKAFSCSS SLSKHKRIHS GDKPYECKEC GKAFSSSSHL IIHIRIHTGE KPYECKECKG AFSESSKLTV HGRTHTGEKP YKCKECKGKAY NCPSSLSIHM RKHTGEKPYE CLECGKAFYL PTSLNTHVKN QSREKPYECK ECGKAFSCPS SFRAHVRDHT GKIQYECKEC GKTFSRSSSL TEHLRTHSGE KPYECKECKG AFISSSHLTV HIRHTHTGEKP YECKKCGKAF IYPSALRIHM RTHHTGEKPYE CKECGKAFRH SSYLTVHARM HTGEKPFECLE ECGKAFSCPS SFRRHVRSHT GEKPYECKEC GKAFVCPAYF RRHVKTHTRE NI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF699
<b>Protein Names</b>	Recommended name: Zinc finger protein 699 Alternative name(s): Hangover homolog
<b>Expression Region</b>	1-642
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.