



# Recombinant Bovine Interferon regulatory factor 1 (IRF1)

<b>Product Code</b>	CSB-EP011814BO-B
<b>Abbreviation</b>	IRF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3SZP0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPITRMRMRP WLEMQINSNQ IPGLIWINKE EMIFQIPWKH AAKHGWDINK DACLFRSWAI HTGRYKAGEK EPDPKTWKAN FRCAMNSLPD IEEVKDQSRN KGSSAVRVYR MLPPLTKSQR KERKSKSSRD ARSKAKKKPY GEYSPDTFSD GLSSSTLPDD HSNYTVRSYM GQDLDIERTL TPALSPCGVS STLPNWSIPV EIVPDSTSDL YNFQVSPMPS TSEAATDEDE EGKLTEDIMK LLEQGTGWQQT SVDGKGYLLN EPGAQPTSVY GEFSCKEEPE VDSPGGYIGL ISSDMKNMDP SLLDSSLTPV RLPSIQAIPC AP
<b>Source</b>	E.coli
<b>Target Names</b>	IRF1
<b>Protein Names</b>	Recommended name: Interferon regulatory factor 1 Short name= IRF-1
<b>Expression Region</b>	1-322
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	IRF1 encodes interferon regulatory factor 1, a member of the interferon regulatory transcription factor (IRF) family. IRF1 serves as an activator of interferons alpha and beta transcription, and in mouse it has been shown to be required for double-stranded RNA induction of these genes. IRF1 also functions as a transcription activator of genes induced by interferons alpha, beta, and gamma. Further, IRF1 has been shown to play roles in regulating apoptosis and tumor-suppressoion.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.