



# Recombinant Human Zinc finger protein 791 (ZNF791)

<b>Product Code</b>	CSB-MP027047HU
<b>Abbreviation</b>	ZNF791
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3KP31
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDSVAFEDVS VSFSQEEWAL LAPSQKKLYR DVMQETFKNL ASIGEKWEDP NVEDQHKNNQG RNLRSHTGER LCEGKEGSQC AENFSPNLSV TKKTAGVKPY ECTICGKAFM RLSSLTRHMR SHTGYELFEK PYKCKECEKA FSYLKSFQRH ERSHTGEKPY KCKQCGKTFI YHQPFRHER THIGEKPYEC KQCGKALSCS SSLRVHERIH TGEKPYECKQ CGKAFSCSSS IRVHERHTG EKPYACKECG KAFISHTSVL THMITHNGDR PYKCKEKGKA FIFPSFLRVH ERIHTGEKPY KCKQCGKAFR CSTSIQIHER IHTGEKPYKC KECGKSFSAR PAFRVHVRVH TGEKPYKCKE CGKAFSRISY FRIHERHTG EKPYECKKCG KTFNYPLDLK IHKRNHTGEK PYECKECAKT FISLENFRRH MITHTGDGPY KCRDCGKVFI FPSALRTHR THTGEKPYEC KQCGKAFSCS SYIRIHKRTH TGEKPYECKE CGKAFIYPTS FQGHMRMHTG EKPYKCKEKG KAFSLHSSFQ RHTRIHNYEK PLECKQCGKA FSVSTSLKKH MRMHNR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF791
<b>Protein Names</b>	Recommended name: Zinc finger protein 791
<b>Expression Region</b>	1-576
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.