



Recombinant Human Ral guanine nucleotide dissociation stimulator-like 3 (RGL3)

Product Code	CSB-EP019626HU-B
Abbreviation	RGL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3MIN7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MERTAGKELA LAPLQDWGEE TEDGAVYSVS LRRQRSQRRS PAEGPGGSQA PSPIANTFLH YRTSKVRVLR AARLERLVGE LVFGDREQDP SFMPAFLATY RTFVPTACLL GFLLPPMPPP PPPGVEIKKT AVQDLSFNKN LRAVVSVLGS WLQDHPQDFR DPPAHS DLGS VRTFLGWAAP GSAEAQKAEK LLEDFLEEAE REQEEPPQV WTGPPRVAQT SDPDSSEACA EEEGLMPQG PQLLDFSVDE VAEQLTLIDL ELFSKVRLYE CLGSVWSQRD RPGAAGASPT VRATVAQFNT VTGCVLGSVL GAPGLAAPQR AQRLEKWIRI AQRCRELRNF SSLRAILSAL QSNPIYRLKR SWGAVSREPL STFRKLSQIF SDENNHLSR EILFQEEATE GSQEEDNTPG SLPSKPPPGP VPYLGTFITD LVMLDTALPD MLEGDLINFE KRRKEWEILA RIQQLRRCQ SYTLSPHPPI LAALHAQNQL TEEQSYRLSR VIEPPAASCP SSPRIRRRIS LTKRLSAKLA REKSSSPSGS PGDPSSPTSS VSPGSPSSP RSRDAPAGSP PASPGPQGPS TKLPLSLDLP SPRPFALPLG SPRIPLPAQQ SSEARVIRVS IDNDHGNLYR SILLTSQDKA PSVRRALQK HNVPQPWACD YQLFQVLP GD RVLLIPDNAN VFYAMSPVAP RDFMLRRKEG TRNTLSVSPS
Source	E.coli
Target Names	RGL3
Protein Names	Recommended name: Ral guanine nucleotide dissociation stimulator-like 3 Short name= RalGDS-like 3
Expression Region	1-710
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.