



Recombinant Human Zinc finger protein 404 (ZNF404)

Product Code	CSB-YP026724HU
Abbreviation	ZNF404
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q494X3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MARVPLTFSD VAIDFSQEEW EYLNSDQRDL YRDVMLENYT NLVSLDFNFT TESNKLSSSEK RNYEVNAYHQ ETWKRNKTFN LMRFIFRTDP QYTIEFGRQQ RPKVGCF SQM IFKKHKSLPL HKRNTREKS YECKEYKKGKGF RKYLHLTEHL RDHTGVIPYE CNECGKAFVV FQHFIRHRKI HTDLKPYECN GCEKAFRFYS QLIQHQIHT GMKPYECKQC GKAFRRHSHL TEHQKIHVGL KPFECKECKE TFRLYRHMCL HQKIHGKVP YKCKECKGAF GHRSSLYQHK KIHSGEKPYK CEQCEKAFVR SYLLVEHQRS HTGKPEHECM ECGKAFSKGS SLLKHKRIHS SEKLYDCKDC GKAFCRGSQL TQHQRHTGE KPHECKECKG TFKLHSYLIQ HQIHTDLKP YECKQCGKAF SRVGD LKTHQ SIHAGEKPYE CKECGKTFRL NSQLIYHQT IHTGLKPYVCK ECKKAFRSIS GLSQHKRIHT GEKPYECKE DKAFNRS DRL TQHETIHTGV KPQKCKECKG AFSHCYQLSQ HQRFHGGERL LM
Source	Yeast
Target Names	ZNF404
Protein Names	Recommended name: Zinc finger protein 404
Expression Region	1-552
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.