



Recombinant Human Zinc finger protein 260 (ZNF260)

Product Code	CSB-BP026620HU
Abbreviation	ZNF260
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3ZCT1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MIGMLESLQH ESDLLQHDQI HTGEKPYECN ECRKTFSLKQ NLVEHKKMHT GEKSHECTEC GKVCSRVSLL TLHLRSHTGK KAYKCNKCGK AFSQKENFLS HQQHHTGEKP YECEKVSQIM PTIIRHQKNH TGTKPYACKE CGKAFNGKAY LTEHEKIHTG EKPFEKNQCG RAFSQKQYLI KHQNIHTGKK PFKCSECGKA FSQKENLIH QRIHTGEKPY ECKGCGKAFI QKSSLIRHQR SHTGEKPYTC KECGKAFSGK SNLTEHEKIH IGEKPYKCNE CGTIFRQKQY LIKHHNIHTG EKPYEKNCCKG KAFSRLSLI VHVRIHTGDK PYECKVCGKA FCQSSSLTVH MRSHTGEKPY GCNECGKAFS QFSTLALHMR IHTGEKPYQC SECGKAFSQK SHHIRHQRIH TH
Source	Baculovirus
Target Names	ZNF260
Protein Names	Recommended name: Zinc finger protein 260 Short name= Zfp-260
Expression Region	1-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.