



Recombinant Bovine Interferon regulatory factor 3 (IRF3)

Product Code	CSB-BP011818BO
Abbreviation	IRF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4JF28
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MGTQKPRILP WLISQLDRGE LEGVAWLGES RTRFRIPWKH GLRQDAQQED FGIFQAWAVA SGAYTPGKDK PDLPTWKRFN RSALNRKEVL RLAEDHSKDS QDPHKIYEFV NSGVRDIPEP DTSQDNGRHN TSQTQEDTLE KLLSDMDLSP GGPSNLTMAS EKPPQFLQSP DSDIPALCPN SGLSENPLKQ LLANEEDWEF EVTAFYRGCQ VFQQTVFCPG GLRLVGSEAG DRMLPGQPIR LPDPATSLTD KSVTDYVQRV LSCLGGGLAL WRAGQWLCAQ RLGHCYVYWA IGEELLPCSG HKPDGEVPKD REGGVFNLGP FITDLITFIE GSRRSPLYTL WFCVGQSWPQ DQPWIKRLVM VKVVPMLLRV LVDIARQGGA SSENTVDLH ISNSDPLSLT PDQYMAQLQD LAEDMDF
Source	Baculovirus
Target Names	IRF3
Protein Names	Recommended name: Interferon regulatory factor 3 Short name= IRF-3
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	IRF3 encodes interferon regulatory factor 3, a member of the interferon regulatory transcription factor (IRF) family. IRF3 is found in an inactive cytoplasmic form that upon serine/threonine phosphorylation forms a complex with CREBBP. This complex translocates to the nucleus and activates the transcription of interferons alpha and beta, as well as other interferon-induced genes.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.