



# Recombinant Bovine Coagulation factor XI (F11), partial

<b>Product Code</b>	CSB-YP007916BO
<b>Abbreviation</b>	F11
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5NTB3
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	EC VTTLFQDACF KGGDITVAFA PNAKHCQIIC THHPRCLLFT FMTESSSEDP TKWYTCILKD SVTETLPMVN MTGAISGYSS KQCLHHISAC SKDMYVDLNM KGMNYNSSLA QSARECQQRQ TDDTHCHFFT FATRHFPSIK DRNTCLLKNT QTGTPTSITK LHEVVSGFSL KSCGLSNLAC IRDIFPRTAF VEDITVMA PDPFVCRSIC THHPSCFFFT FLSEEWPTAS ERNLCLLKTS SSSLPSARFR KNRAFSGFSL QHCQHSVPVF CHSSFYRNTD FLGEELDIVD ADSHEACQKT CTNSIRCQFF TYSPSQESCN GGKGGKCYLKL SANGSPTKIL HGTGSISGYT LRLCKMDNVC TTKIKTR
<b>Source</b>	Yeast
<b>Target Names</b>	F11
<b>Protein Names</b>	Recommended name: Coagulation factor XI Short name= FXI EC= 3.4.21.27Alternative name(s): Plasma thromboplastin antecedent Short name= PTACleaved into the following 2 chains: 1. Coagulation factor XIa heavy chain 2.
<b>Expression Region</b>	19-387
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.