



Recombinant *Oryza sativa* subsp. japonica Mitogen-activated protein kinase 15 (MPK15)

Product Code	CSB-MP013454OFG
Abbreviation	MPK15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53N72
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDFFTEYGEG NRYKIEEVIG KGSYGVVCSA LDHTTGEKVA IKKINDIFEH VSDATRILRE IKLLRLLRHP DIVEIKHILL PPSRREFKDI YVVFELMESD LHQVIKANDD LTPEHYQFFL YQLLRGLKYI HTANVFHRDL KPKNILANAD CKLKICDFGL ARVAFSDTPT AIFWTDYVAT RWYRAPELPG SFFSKYTPAI DIWSIGCIFA ELLTGKPLFP GKNVVHQLDI ITDLLGTPST EAISRIRNEK ARRYLSSMRR KKPIPFQKF PNADPLALRL LERMLSFEPEK DRPNAAEEALA DPYFRNIANV DREPSAQPVV KLEFEFERRR ITKEDIRELI YRDILEYHPN MLREYLEGTE SAGFMYPYSAV DHFKKQFAYL EEHYAKGSTA APPERQHNSL PRPSVLYSDD RPQNTANIAE DLSKCVLGDN TQKMHQGSAS VCANRVPQGG AARPGKVVGS ALRYGNCSTS TAEQYEHRRR DRNPALATNT VSPRGSYP
Source	Mammalian cell
Target Names	MPK15
Protein Names	Recommended name: Mitogen-activated protein kinase 15 Short name= MAP kinase 15 EC= 2.7.11.24
Expression Region	1-498
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.