



# Recombinant Pig Transcription factor A, mitochondrial (TFAM)

<b>Product Code</b>	CSB-EP023413PI-B
<b>Abbreviation</b>	TFAM
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5D144
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SSTLSGFPPKPMTSYVRFSEQLPIFKAQNPDAKNSELIKKIAELWRELPDSEK KIYEDA YRADWQVYKEEVNRIQEQLTPSQMVSLEKEIMQKRLKKKALIKKRELTMLGKP KRPRSAY NIFIAERFQEAKDGPSQVKLKTINENWKNLSSSQKQVYIQLAEDDKVRYYNEM KSWEEQM VEVGRNDLIRRSMKHSAKKDTEEC
<b>Source</b>	E.coli
<b>Target Names</b>	TFAM
<b>Protein Names</b>	Recommended name: Transcription factor A, mitochondrial Short name= mtTFA
<b>Expression Region</b>	43-246
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a mitochondrial transcription factor that is a key activator of mitochondrial transcription as well as a participant in mitochondrial genome replication. Studies in mice have demonstrated that this gene product is required to regulate the mitochondrial genome copy number and is essential for embryonic development. A mouse model for Kearns-Sayre syndrome was produced when expression of this gene was eliminated by targeted disruption in heart and muscle cells.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.