



Recombinant *Saccharomyces cerevisiae* Putative lipase YOR059C (YOR059C)

Product Code	CSB-EP600105SVG-B
Abbreviation	YOR059C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q08448
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MTSDKHLFVLIHGLWGN YTHMESMRTILSTTLKKEDVNDDMIYFLPKQNAMFK TFDGIEI IGYRTLIEVCEFIRDYKDGKITKLSVMGYSQGGGLVARFMIGKMLTEFKELFEDIE PQLFI TMATPHLGVEFYNPTGIA YKSALYSALRTL GSTILGKSGREMFIA NSNNILVKL SQGEY LEALSLFKWRIAFANVKNDR TVAFYTA FITDCDPFIDFNKLYTFEEKIPGSGY KGILP KIVDLNALNVNSHAPTKPTKYKKWGR TILILVATFLILPIALVMNGLGTAYSIV TCK YRKMLSN GILHNEVRGKLGLTEQLKGYVTDAYGSIINSALDMDANYEASNSNL VNEEELP WKEFIQKYTTINDGVWKSFKKLPFDENRKVILRNLNKLKWIRVPIYIKAVNAHG VIVAR RGM DENTAATGIACIEFTAQLLAYLMHKS N
Source	E.coli
Target Names	YOR059C
Protein Names	Recommended name: Putative lipase YOR059C EC= 3.1.-.-
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.