



Recombinant Schizosaccharomyces pombe 5-azacytidine resistance protein azr1 (azr1)

Product Code	CSB-BP600213SXV
Abbreviation	azr1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09189
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MNVTRLYFRV AGTKQFARYV HKYAAYSSTS FQKKKSHFPS PATLDHPDAG EDAFINLRNE NYILNAVFDG VGGWANVGID PSIFSWGLVR EIKKVFNNSD EFQPSPLTLL SKAYAALKKS NTVEAGSSTA CLTLFNCNG KLHSLNLGDS GFLILRNGAI HYASPAQLQ FNMPYQLAIY PRNYRSAENI GPKMGQATVH DLKDNDLVIL ATDGIFDNIE EKSILDIAGV VDFSSLSNVQ KCLDDLAMRI CRQAVLNSLD TKWESPFAKT AKSFGFKFQG GKVDDTTITC LLIKSKNYE
Source	Baculovirus
Target Names	azr1
Protein Names	Recommended name: 5-azacytidine resistance protein azr1
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.