



Recombinant Probable dynactin subunit 2 (dnc-2)

Product Code	CSB-MP600218CXY
Abbreviation	dnc-2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q09248
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	<p>MSSIGKEDIY ESDGAASQVQ KSEPENDNPT HPDIELISVD VDEALKKYKN RVLNLTSSDF SDSIAKKRRH AFGNNQYVLE VTGAGFSGTE TAAEKLNRIL YEVA DLNEQI RADENLKTDL LNAEVLENLE KEVKTQLQVAQ SNGKTARVEH DVLPNVRTD SKVATLENRL RRIEQVIGSS VIPAPVLDI IEDLKLRCET LNHSYVSGLE QRLNMMMLTKL EKIDETRANN DIDENLDKKV DEILELMQKW DVACSSLPSN VNKVKALNRL HEQVQHFAGG LSHLKTIREK LEKEVAQGRE AIIYESQGK QEIGSVVEKL KLLLEAKVADL K</p>
Source	Mammalian cell
Target Names	dnc-2
Protein Names	Recommended name: Probable dynactin subunit 2
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>