



Recombinant Uncharacterized protein C32A3.2 (C32A3.2)

Product Code	CSB-YP600220CXY
Abbreviation	C32A3.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09261
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	SSAASGDGN DGKGAKDAID EDLLNAIEGV ANNIHPQNGS EKKSLKNTLI NRLVANEKAS FAAAAASASS EMLDDQALIG LLADVAGDAK VEKKLPPKSA QLRQEKRGVLV LLRKEIFYQA VQSGFTTEEA RVKSETIVNE AQIKLQEQRK ALLNDVREKV EQEEVEETER SEKDQKLFTM ALEFMEKIYK DDLISSAVRK PVKVDNDAIK LFNQKPLGIW KKGEKYEDFS LGFWKQWDER AARISNGSFG PTNSFEEQIE WTTKGKQWEY PIDNEFKMGD ESNVSFIDHV FLERHLPSLG IPKSGPIAHF MHLVCVGLSK NPYMTAAKKR EHLKWFADYF NTEKQKL VHK LHEQE QIAAQ NAL
Source	Yeast
Target Names	mrps-31
Protein Names	Recommended name: Uncharacterized protein C32A3.2
Expression Region	22-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.