



Recombinant Splicing factor 3B subunit 4 (sap-49)

Product Code	CSB-BP600239CXY
Abbreviation	sap-49
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09442
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSAGPIVERN QDATIYVGGL DEKVSESILW ELMVQAGPVV SVNMPKDRVT ANHQGFGFVE FMGEEDADYA IKILNMIKLY GKPIKVNKAS AHEKNMDVGA NIFVGNLDPE VDEKLLYDTF SAFGVILQVP KIMRDVDSGT SKGFAFINFA SFEASDTALE AMNGQFLCNR AITVSYAFKR DSKGERHGTA AERMLAAQNP LFPKDRPHQV FSDVPLGVPA NTPLAMPGVH AAIAAHATGR PGYQPPPLMG MAQSGYQGQY PPVPPPPSV TPMPPMPPT PGMTPRPPPP PSSGMWPPPP PPPGRTPGP PGMPGMPPPP PPSRFGPPGM GGMPPPPPG MRYPPGMPPP PPRYPSAGP GMYPPPPSR PPAPPSGHGM IPPPPPS
Source	Baculovirus
Target Names	sap-49
Protein Names	Recommended name: Splicing factor 3B subunit 4 Alternative name(s): Spliceosome-associated protein 49
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.