



Recombinant Probable dimethyladenosine transferase (E02H1.1)

Product Code	CSB-BP600247CXY
Abbreviation	E02H1.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09522
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MGKTSKVKKT KAGSSTGNVQ SLPFNTDKGQ HILKNPGVVN AIVEKSALKA TDTVLEVGPG TGNLTVKMLE VAKTVIACEI DPRMIAEVKK RVMGTPLQNK LQVNGGDVMK MEWPFVDVCV ANLPYQISSP FVQKLLLHRP LPRYAVLMFQ KEFADRLVAR PGDKDYSRLS VNVQLLAKVE MLMKVKRTEF RPPPKVDSAV VRIAPKNPPP PVNFVEWEGL LRLCFMRKNK TLMAIFRLSN VIEVIEDNFR KVCSFKNKPI PKDLNMKKVI EETLTASGYG ESRARKMRVE DFLALLLAFN KADIHFLS
Source	Baculovirus
Target Names	E02H1.1
Protein Names	Recommended name: Probable dimethyladenosine transferase EC= 2.1.1.183 Alternative name(s): Probable 18S rRNA (adenine(1779)-N(6)/adenine(1780)-N(6))-dimethyltransferase Probable 18S rRNA dimethylase Probable S-adenosylmethionine-6-
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.