



Recombinant ATP synthase subunit delta, mitochondrial (F58F12.1)

Product Code	CSB-EP600249CXY
Abbreviation	F58F12.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09544
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	AA APAANANPEE LRLTFASPDV AVFSNAVVKQ VDVPTLAGMV GVLANHVPIT GVLKPGVSVV TTNEGTVQRL FVSSGTL SVN IDGSCQVLAE EVLKVEEIDE SAARAELDAA QRASGEGSEV ARAEAQIRAE VAEALIKAAT NQQ
Source	E.coli
Target Names	F58F12.1
Protein Names	Recommended name: ATP synthase subunit delta, mitochondrial Alternative name(s): F-ATPase delta subunit
Expression Region	19-163
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.