



Recombinant Proteasome subunit alpha type-3 (pas-7)

Product Code	CSB-BP600252CXY
Abbreviation	pas-7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q09583
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	SSIGTGYDL AASTFSPDGR IFQVEYAQKA VDNAGTMIAI RGKNGVVVVA DKLISSKLYT DNANPRMFNV NDNVGVAVAG NYPDGFALKN YAYGEAMKWL KDYREMPIQ NIANVAEYI HIHTLGISRP FGAGAFFMSW NKQTGGRLFL VEPSGLNYEY KAWAVGKHRQ AAKAEIEKLE IEELDVNQLV KEAARIIMVV RDENKDKNVQ IEMGWVGEQT NGKYEEVPSE VVTAEEWAI AKLDEDDMDD
Source	Baculovirus
Target Names	pas-7
Protein Names	Recommended name: Proteasome subunit alpha type-3 EC= 3.4.25.1 Alternative name(s): Proteasome subunit alpha 7
Expression Region	2-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.